

GerMed Periotomes Kit Set

*GerDent*USA
www.GerDentUSA.com

If we don't have the instrument
you're looking for,
we will **MANUFACTURE** it for you!

The GerMed Periotomes Kit set has 8 unique periotomes for precisely cutting PDL before tooth extractions. Each instrument features unique blades for targeting tooth structures with minimal chances of complications and damage to the surrounding tissues.



SKU: GD50-1445



☎ 516 593 7100 - 516 807 0884
📍 61-21 Springfield Blvd, Suite 396 Oakland
Gardens, NY 11364 United States

GerMed Periotomes Kit Set

SKU: GD50-1436

GerMed Periotome #1 PT4

Used for atraumatic tooth extraction by carefully severing the periodontal ligament.

- Sharp, narrow blade
- PT4 configuration
- German stainless steel
- Ergonomic handle



SKU: GD50-1437

GerMed Periotome #2 PT5

Specially designed to cut periodontal ligaments during tooth extraction procedures.

- PT5 Configuration
- Sharp Tip
- Lightweight



SKU: GD50-1438

GerMed Periotome #3

Used for atraumatic tooth extraction by severing the periodontal ligament with precision.

- Sharp, thin blade
- Sharp tip
- German stainless steel
- Ergonomic handle



SKU: GD50-1439

GerMed Periotome #4 PT3

Specialized dental instrument designed to section or divide periodontal ligaments during tooth extraction procedures.

- PT3 Configuration
- Sharp Tip
- German Stainless Steel




SKU: GD50-1440

GerMed Periotome #5 Small Gouge

specialized dental instrument designed to section or divide periodontal ligaments during tooth extraction procedures.

- Small Gouge Shape
- Sharp Tip
- Surgical-grade stainless steel



SKU: GD50-1442

GerMed Periotome #6

specialized dental instrument used to cut periodontal ligaments during tooth extraction procedures.

- Precision Blade
- Sharp Tip
- Ergonomic Design
- German stainless steel



SKU: GD50-1443

GerMed Periotome #7

specialized dental instrument used to section or divide periodontal ligaments during tooth extraction procedures.

- Precision Blade
- Sharp Tip
- Ergonomic Design



516 593 7100
516 807 0884



www.GerDentUSA.com